

Recombinant Human MARK4 Protein, N-His

Catalog No: #AP97744

Package Size: #AP97744-1 50ug #AP97744-2 100ug #AP97744-3 1mg

Description

Product Name	Recombinant Human MARK4 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	MAP/microtubule affinity-regulating kinase 4, 2.7.11.1, MAP/microtubule affinity-regulating kinase-like 1, MARK4, KIAA1860, MARKL1
Accession No.	Q96L34
Target Species	Homo sapiens (Human)
Calculated MW	38.41 kDa
Target Length	Gly57-Lys370
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.