

Recombinant Human HDGFL2 Protein, N-His-SUMO & C-Strep

Catalog No: #AP97748

Package Size: #AP97748-1 50ug #AP97748-2 100ug #AP97748-3 1mg

Description

Product Name	Recombinant Human HDGFL2 Protein, N-His-SUMO & C-Strep
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Hepatoma-derived growth factor-related protein 2, HDGF-related protein 2, HRP-2, Hepatoma-derived growth factor 2, HDGF-2, HDGFL2, HDGF2, HDGFRP2, HRP2, UNQ785/PRO1604
Accession No.	Q7Z4V5
Target Species	Homo sapiens (Human)
Calculated MW	23.62 kDa
Target Length	Met1-Asn87
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.