

Recombinant Human USP33 Protein, N-His-SUMO & C-Strep

Catalog No: #AP97754

Package Size: #AP97754-1 50ug #AP97754-2 100ug #AP97754-3 1mg

Description

Product Name	Recombinant Human USP33 Protein, N-His-SUMO & C-Strep
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Ubiquitin carboxyl-terminal hydrolase 33, 3.4.19.12, Deubiquitinating enzyme 33, Ubiquitin thioesterase 33, Ubiquitin-specific-processing protease 33, VHL-interacting deubiquitinating enzyme 1, hVDU1, USP33, KIAA1097, VDU1
Accession No.	Q8TEY7
Target Species	Homo sapiens (Human)
Calculated MW	23.23 kDa
Target Length	Arg36-Glu120
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.