

Recombinant Human POFUT1 Protein, N-GST & C-His

Catalog No: #AP97771

Package Size: #AP97771-1 50ug #AP97771-2 100ug #AP97771-3 1mg

Description

Product Name	Recombinant Human POFUT1 Protein, N-GST & C-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	GDP-fucose protein O-fucosyltransferase 1, 2.4.1.221, Peptide-O-fucosyltransferase 1, O-FucT-1, POFUT1, FUT12, KIAA0180
Accession No.	Q9H488
Target Species	Homo sapiens (Human)
Calculated MW	39.02 kDa
Target Length	Leu284-Asp379
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.