

Recombinant Human IRGM Protein, N-His

Catalog No: #AP97773

Package Size: #AP97773-1 50ug #AP97773-2 100ug #AP97773-3 1mg

Description

Product Name	Recombinant Human IRGM Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Immunity-related GTPase family M protein, 3.6.5.-, Immunity-related GTPase family M protein 1, Interferon-inducible protein 1, LPS-stimulated RAW 264.7 macrophage protein 47 homolog, LRG-47, IRGM, IFI1, IRGM1, LRG47
Accession No.	A1A4Y4
Target Species	Homo sapiens (Human)
Calculated MW	22.31 kDa
Target Length	Met1-Tyr181
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.