

Recombinant Human TARSS3 Protein, N-His

Catalog No: #AP97776

Package Size: #AP97776-1 50ug #AP97776-2 100ug #AP97776-3 1mg

Description

Product Name	Recombinant Human TARSS3 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Threonine--tRNA ligase 2, cytoplasmic, 6.1.1.3, Threonyl-tRNA synthetase, ThrRS, Threonyl-tRNA synthetase protein 3, TARSS3, TARSL2
Accession No.	A2RTX5
Target Species	Homo sapiens (Human)
Calculated MW	28.01 kDa
Target Length	Thr158-Asp383
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.