

Recombinant Human BDP1 Protein, N-GST & C-His

Catalog No: #AP97779

Package Size: #AP97779-1 50ug #AP97779-2 100ug #AP97779-3 1mg

Description

Product Name	Recombinant Human BDP1 Protein, N-GST & C-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Transcription factor TFIIIB component B" homolog, Transcription factor IIIB 150, TFIIIB150, Transcription factor-like nuclear regulator, BDP1, KIAA1241, KIAA1689, TFNR
Accession No.	A6H8Y1
Target Species	Homo sapiens (Human)
Calculated MW	39.82 kDa
Target Length	Ser300-Ser396
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.