

Recombinant Human SLC9A9 Protein, N-His-KSI

Catalog No: #AP97784

Package Size: #AP97784-1 50ug #AP97784-2 100ug #AP97784-3 1mg

Description

Product Name	Recombinant Human SLC9A9 Protein, N-His-KSI
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Sodium/hydrogen exchanger 9, Na(+)/H(+) exchanger 9, NHE-9, Solute carrier family 9 member 9, SLC9A9, NHE9, Nbla00118
Accession No.	Q8IVB4
Target Species	Homo sapiens (Human)
Calculated MW	40.55 kDa
Target Length	[Arg67-Lys126] & [Leu487-Asn645]
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.