

Recombinant Human KDM4E Protein, N-His

Catalog No: #AP97786

Package Size: #AP97786-1 50ug #AP97786-2 100ug #AP97786-3 1mg

Description

Product Name	Recombinant Human KDM4E Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Lysine-specific demethylase 4E, 1.14.11.66, KDM4D-like protein, Lysine-specific demethylase 4D-like, [histone H3]-trimethyl-L-lysine(9) demethylase 4E, KDM4E, KDM4DL
Accession No.	B2RXH2
Target Species	Homo sapiens (Human)
Calculated MW	43.61 kDa
Target Length	Thr14-Arg372
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.