

Recombinant Human Lambda chain Protein, N-His

Catalog No: #AP97790

Package Size: #AP97790-1 50ug #AP97790-2 100ug #AP97790-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Lambda chain Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Immunoglobulin lambda constant 1, Ig lambda chain C region MGC, Ig lambda-1 chain C region, IGLC1
Accession No.	P0CG04
Target Species	Homo sapiens (Human)
Calculated MW	12.50 kDa
Target Length	Gly1-Thr95
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.