

Recombinant Human ECRG4 Protein, N-His-KSI

Catalog No: #AP97791

Package Size: #AP97791-1 50ug #AP97791-2 100ug #AP97791-3 1mg

Description

Product Name	Recombinant Human ECRG4 Protein, N-His-KSI
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Augurin, Esophageal cancer-related gene 4 protein, ECRG4, ECRG4, C2orf40, UNQ761/PRO1508
Accession No.	Q9H1Z8
Target Species	Homo sapiens (Human)
Calculated MW	29.97 kDa
Target Length	Asn32-Tyr148
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.