

Recombinant Human GPBAR1 Protein, N-His-KSI

Catalog No: #AP97796

Package Size: #AP97796-1 50ug #AP97796-2 100ug #AP97796-3 1mg

Description

Product Name	Recombinant Human GPBAR1 Protein, N-His-KSI
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	G-protein coupled bile acid receptor 1, G-protein coupled receptor GPCR19, hGPCR19, Membrane-type receptor for bile acids, M-BAR, hBG37, BG37, GPBAR1, TGR5
Accession No.	Q8TDU6
Target Species	Homo sapiens (Human)
Calculated MW	28.29 kDa
Target Length	[Thr2-Gly19] & [Pro72-Cys85] & [Asn147-Tyr165] & [Ala250-Thr261] & [Gly283-Asn330]
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.