

Recombinant Mouse CD3D Protein, N-GST & C-His

Catalog No: #AP97803

Package Size: #AP97803-1 50ug #AP97803-2 100ug #AP97803-3 1mg

Description

Product Name	Recombinant Mouse CD3D Protein, N-GST & C-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	T-cell surface glycoprotein CD3 delta chain, T-cell receptor T3 delta chain, CD3d, Cd3d, T3d
Accession No.	P04235
Target Species	Mus musculus (Mouse)
Calculated MW	37.54 kDa
Target Length	Phe22-Ala105
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.