

Recombinant Mouse CD54/ICAM1 Protein, N-His

Catalog No: #AP97804

Package Size: #AP97804-1 50ug #AP97804-2 100ug #AP97804-3 1mg

Description

Product Name	Recombinant Mouse CD54/ICAM1 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Intercellular adhesion molecule 1, ICAM-1, MALA-2, MyD10, CD54, Icam1, Icam-1
Accession No.	P13597
Target Species	Mus musculus (Mouse)
Calculated MW	24.40 kDa
Target Length	Leu120-Gln321
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.