

Recombinant Mouse CD126/IL6R/IL-6RA Protein, N-His

Catalog No: #AP97806

Package Size: #AP97806-1 50ug #AP97806-2 100ug #AP97806-3 1mg

Description

Product Name	Recombinant Mouse CD126/IL6R/IL-6RA Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Interleukin-6 receptor subunit alpha, IL-6 receptor subunit alpha, IL-6R subunit alpha, IL-6R-alpha, IL-6RA, IL-6R 1, CD126, Soluble interleukin-6 receptor subunit alpha, sIL6R, Il6ra, Il6r
Accession No.	P22272
Target Species	Mus musculus (Mouse)
Calculated MW	33.26 kDa
Target Length	Gly33-Thr310
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.