

Recombinant Mouse CD224/GGT1 Protein, N-His

Catalog No: #AP97812

Package Size: #AP97812-1 50ug #AP97812-2 100ug #AP97812-3 1mg

Description

Product Name	Recombinant Mouse CD224/GGT1 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Glutathione hydrolase 1 proenzyme, 3.4.19.13, Gamma-glutamyltransferase 1, Gamma-glutamyltranspeptidase 1, GGT 1, 2.3.2.2, Leukotriene-C4 hydrolase, 3.4.19.14, CD224, Glutathione hydrolase 1 heavy chain, Glutathione hydrolase 1 light chain, Ggt1, Ggt, Ggtp
Accession No.	Q60928
Target Species	Mus musculus (Mouse)
Calculated MW	59.61 kDa
Target Length	Arg39-Tyr568
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.