

Recombinant Mouse CD11c/ITGAX Protein, N-His

Catalog No: #AP97813

Package Size: #AP97813-1 50ug #AP97813-2 100ug #AP97813-3 1mg

Description

Product Name	Recombinant Mouse CD11c/ITGAX Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Integrin alpha-X, CD11 antigen-like family member C, Leukocyte adhesion glycoprotein p150,95 alpha chain, Leukocyte adhesion receptor p150,95, CD11c, Itgax
Accession No.	Q9QXH4
Target Species	Mus musculus (Mouse)
Calculated MW	23.91 kDa
Target Length	Asp150-Thr342
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.