

Recombinant Mouse CD274/PD-L1/B7-H1 Protein, N-His

Catalog No: #AP97824

Package Size: #AP97824-1 50ug #AP97824-2 100ug #AP97824-3 1mg

Description

Product Name	Recombinant Mouse CD274/PD-L1/B7-H1 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Programmed cell death 1 ligand 1, PD-L1, PDCD1 ligand 1, Programmed death ligand 1, B7 homolog 1, B7-H1, CD274, Cd274, B7h1, Pdcd1l1, Pdcd1lg1, Pdl1
Accession No.	Q9EP73
Target Species	Mus musculus (Mouse)
Calculated MW	27.19 kDa
Target Length	Phe19-His239
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.