

Recombinant Mouse MASP2 Protein, N-His

Catalog No: #AP97825

Package Size: #AP97825-1 50ug #AP97825-2 100ug #AP97825-3 1mg

Description

Product Name	Recombinant Mouse MASP2 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Mannan-binding lectin serine protease 2, 3.4.21.104, MBL-associated serine protease 2, Mannose-binding protein-associated serine protease 2, MASP-2, Mannan-binding lectin serine protease 2 A chain, Mannan-binding lectin serine protease 2 B chain, Masp2
Accession No.	Q91WP0
Target Species	Mus musculus (Mouse)
Calculated MW	47.60 kDa
Target Length	Glu24-Glu430
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.