

Recombinant Mouse AIM2 Protein, N-His

Catalog No: #AP97827

Package Size: #AP97827-1 50ug #AP97827-2 100ug #AP97827-3 1mg

Description

Product Name	Recombinant Mouse AIM2 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Interferon-inducible protein AIM2, Interferon-inducible protein 210, Ifi-210, Interferon-inducible protein p210, Aim2, Gm1313, Ifi210
Accession No.	Q91VJ1
Target Species	Mus musculus (Mouse)
Calculated MW	23.93 kDa
Target Length	Glu153-Lys342
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.