

Recombinant Mouse TMPRSS2 Protein, N-His

Catalog No: #AP97828

Package Size: #AP97828-1 50ug #AP97828-2 100ug #AP97828-3 1mg

Description

Product Name	Recombinant Mouse TMPRSS2 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Transmembrane protease serine 2, 3.4.21.122, Epitheliasin, Plasmic transmembrane protein X, Transmembrane protease serine 2 non-catalytic chain, Transmembrane protease serine 2 catalytic chain, Tmprss2
Accession No.	Q9JIQ8
Target Species	Mus musculus (Mouse)
Calculated MW	43.84 kDa
Target Length	Thr114-Ser490
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.