

Recombinant Mouse H6PD Protein, N-His

Catalog No: #AP97835

Package Size: #AP97835-1 50ug #AP97835-2 100ug #AP97835-3 1mg

Description

Product Name	Recombinant Mouse H6PD Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	GDH/6PGL endoplasmic bifunctional protein, Hexose-6-phosphate dehydrogenase, Glucose 1-dehydrogenase, 1.1.1.47, Glucose-6-phosphate dehydrogenase, 1.1.1.363, 6-phosphogluconolactonase, 6PGL, 3.1.1.31, H6pd
Accession No.	Q8CFX1
Target Species	Mus musculus (Mouse)
Calculated MW	28.68 kDa
Target Length	Ser551-Gly789
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.