

Recombinant Mouse PLAT/tPA Protein, N-His

Catalog No: #AP97836

Package Size: #AP97836-1 50ug #AP97836-2 100ug #AP97836-3 1mg

Description

Product Name	Recombinant Mouse PLAT/tPA Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Tissue-type plasminogen activator, t-PA, t-plasminogen activator, tPA, 3.4.21.68, Tissue-type plasminogen activator chain A, Tissue-type plasminogen activator chain B, Plat
Accession No.	P11214
Target Species	Mus musculus (Mouse)
Calculated MW	29.57 kDa
Target Length	Asp316-Gln559
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.