

Recombinant Mouse GNRH1 Protein, N-GST & C-His

Catalog No: #AP97840

Package Size: #AP97840-1 50ug #AP97840-2 100ug #AP97840-3 1mg

Description

Product Name	Recombinant Mouse GNRH1 Protein, N-GST & C-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Progonadoliberin-1, Progonadoliberin I, Gonadoliberin-1, Gonadoliberin I, Gonadotropin-releasing hormone I, GnRH-I, Luliberin I, Luteinizing hormone-releasing hormone I, LH-RH I, Prolactin release-inhibiting factor 1, Prolactin release-inhibiting factor I, GnRH1, GnRH
Accession No.	P13562
Target Species	Mus musculus (Mouse)
Calculated MW	36.33 kDa
Target Length	Gln22-Met90
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.