

Recombinant Mouse APOH Protein, N-GST & C-His

Catalog No: #AP97846

Package Size: #AP97846-1 50ug #AP97846-2 100ug #AP97846-3 1mg

Description

Product Name	Recombinant Mouse APOH Protein, N-GST & C-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Beta-2-glycoprotein 1, APC inhibitor, Activated protein C-binding protein, Apolipoprotein H, Apo-H, Beta-2-glycoprotein I, B2GPI, Beta(2)GPI, Apoh, B2gp1
Accession No.	Q01339
Target Species	Mus musculus (Mouse)
Calculated MW	37.74 kDa
Target Length	Ser263-Cys345
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.