

Recombinant Mouse F11 Protein, N-His

Catalog No: #AP97851

Package Size: #AP97851-1 50ug #AP97851-2 100ug #AP97851-3 1mg

Description

Product Name	Recombinant Mouse F11 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Coagulation factor XI, FXI, 3.4.21.27, Plasma thromboplastin antecedent, PTA, Coagulation factor XIa heavy chain, Coagulation factor XIa light chain, F11
Accession No.	Q91Y47
Target Species	Mus musculus (Mouse)
Calculated MW	27.91 kDa
Target Length	Ala395-Val624
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.