

Recombinant Mouse S100A9 Protein, N-His

Catalog No: #AP97859

Package Size: #AP97859-1 50ug #AP97859-2 100ug #AP97859-3 1mg

Description

Product Name	Recombinant Mouse S100A9 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Protein S100-A9, Calgranulin-B, Leukocyte L1 complex heavy chain, Migration inhibitory factor-related protein 14, MRP-14, p14, S100 calcium-binding protein A9, S100a9, Cagb, Mrp14
Accession No.	P31725
Target Species	Mus musculus (Mouse)
Calculated MW	15.21 kDa
Target Length	Ala2-Lys113
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.