

Recombinant Mouse GPX1 Protein, N-His

Catalog No: #AP97863

Package Size: #AP97863-1 50ug #AP97863-2 100ug #AP97863-3 1mg

Description

Product Name	Recombinant Mouse GPX1 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Glutathione peroxidase 1, GPx-1, GSHPx-1, 1.11.1.9, Cellular glutathione peroxidase, Phospholipid-hydroperoxide glutathione peroxidase GPX1, 1.11.1.12, Selenium-dependent glutathione peroxidase 1, Gpx1
Accession No.	P11352
Target Species	Mus musculus (Mouse)
Calculated MW	19.81 kDa
Target Length	Gly48-Ser201
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.