

Recombinant Mouse Cyp21 Protein, N-His

Catalog No: #AP97866

Package Size: #AP97866-1 50ug #AP97866-2 100ug #AP97866-3 1mg

Description

Product Name	Recombinant Mouse Cyp21 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Steroid 21-hydroxylase, 1.14.14.16, 21-OHase, Cytochrome P-450c21, Cytochrome P450 21, Cytochrome P450 XXI, Cytochrome P450-C21, Cyp21, Cyp21a-1, Cyp21a1
Accession No.	P03940
Target Species	Mus musculus (Mouse)
Calculated MW	52.98 kDa
Target Length	Leu30-Arg476
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.