

Recombinant Mouse C4B Protein, N-His

Catalog No: #AP97868

Package Size: #AP97868-1 50ug #AP97868-2 100ug #AP97868-3 1mg

Description

Product Name	Recombinant Mouse C4B Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Complement C4-B, Complement C4 beta chain, Complement C4 alpha chain, C4a anaphylatoxin, Complement C4 gamma chain, C4b, C4
Accession No.	P01029
Target Species	Mus musculus (Mouse)
Calculated MW	38.10 kDa
Target Length	Glu988-Thr1317
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.