

Recombinant Mouse SCGB1A1/Uteroglobulin Protein, N-GST & C-His

Catalog No: #AP97869

Package Size: #AP97869-1 50ug #AP97869-2 100ug #AP97869-3 1mg

Description

Product Name	Recombinant Mouse SCGB1A1/Uteroglobulin Protein, N-GST & C-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Uteroglobulin, Club cell 17 kDa protein, Club cell phospholipid-binding protein, CCPBP, Club cells 10 kDa secretory protein, CC10, PCB-binding protein, Secretoglobin family 1A member 1, Scgb1a1, Cc10, Ugb, Utg
Accession No.	Q06318
Target Species	Mus musculus (Mouse)
Calculated MW	38.57 kDa
Target Length	Lys2-Phe96
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.