

Recombinant Mouse ALOX12 Protein, N-His

Catalog No: #AP97873

Package Size: #AP97873-1 50ug #AP97873-2 100ug #AP97873-3 1mg

Description

Product Name	Recombinant Mouse ALOX12 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Polyunsaturated fatty acid lipoyxygenase ALOX12, 1.13.11.-, Arachidonate (12S)-lipoyxygenase, 12S-LOX, 12S-lipoyxygenase, 1.13.11.31, Arachidonate (15S)-lipoyxygenase, 1.13.11.33, Linoleate (13S)-lipoyxygenase, Lipoxin synthase 12-LO, 3.3.2.-, Platelet-type lipoyxygenase 12, P-12LO, Alox12, Alox12p
Accession No.	P39655
Target Species	Mus musculus (Mouse)
Calculated MW	21.35 kDa
Target Length	Met1-His167
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.