

Recombinant Mouse ENPP1 Protein, N-His

Catalog No: #AP97876

Package Size: #AP97876-1 50ug #AP97876-2 100ug #AP97876-3 1mg

Description

Product Name	Recombinant Mouse ENPP1 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Ectonucleotide pyrophosphatase/phosphodiesterase family member 1, E-NPP 1, Lymphocyte antigen 41, Ly-41, Phosphodiesterase I/nucleotide pyrophosphatase 1, Plasma-cell membrane glycoprotein PC-1, Alkaline phosphodiesterase I, 3.1.4.1, Nucleotide pyrophosphatase, NPPase, 3.6.1.9, Nucleotide diphosphatase, Ectonucleotide pyrophosphatase/phosphodiesterase family member 1, secreted form, Enpp1, Npps, Pc1, Pdnpp1
Accession No.	P06802
Target Species	Mus musculus (Mouse)
Calculated MW	36.29 kDa
Target Length	Lys617-Asp906
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.