

Recombinant Mouse CASP1/Caspase-1 Protein, N-His

Catalog No: #AP97881

Package Size: #AP97881-1 50ug #AP97881-2 100ug #AP97881-3 1mg

Description

Product Name	Recombinant Mouse CASP1/Caspase-1 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Caspase-1, CASP-1, 3.4.22.36, Interleukin-1 beta convertase, IL-1BC, Interleukin-1 beta-converting enzyme, ICE, IL-1 beta-converting enzyme, p45, Caspase-1 subunit p20, Caspase-1 subunit p10, Casp1, Il1bc
Accession No.	P29452
Target Species	Mus musculus (Mouse)
Calculated MW	32.97 kDa
Target Length	Leu137-His402
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.