

Recombinant Mouse CDKN2A/p16INK4a Protein, N-His

Catalog No: #AP97888

Package Size: #AP97888-1 50ug #AP97888-2 100ug #AP97888-3 1mg

Description

Product Name	Recombinant Mouse CDKN2A/p16INK4a Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Cyclin-dependent kinase inhibitor 2A, Cyclin-dependent kinase 4 inhibitor A, CDK4I, p16-INK4a, p16-INK4, Cdkn2a, P16ink4a
Accession No.	P51480
Target Species	Mus musculus (Mouse)
Calculated MW	16.64 kDa
Target Length	Met1-Ser135
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.