

Recombinant Mouse BCAR1 Protein, N-GST & C-His

Catalog No: #AP97897

Package Size: #AP97897-1 50ug #AP97897-2 100ug #AP97897-3 1mg

Description

Product Name	Recombinant Mouse BCAR1 Protein, N-GST & C-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Breast cancer anti-estrogen resistance protein 1, CRK-associated substrate, p130cas, Bcar1, Cas, Crkas
Accession No.	Q61140
Target Species	Mus musculus (Mouse)
Calculated MW	35.38 kDa
Target Length	Thr2-Met66
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.