

Recombinant Mouse MECOM Protein, N-His

Catalog No: #AP97903

Package Size: #AP97903-1 50ug #AP97903-2 100ug #AP97903-3 1mg

Description

Product Name	Recombinant Mouse MECOM Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Histone-lysine N-methyltransferase MECOM, 2.1.1.367, Ecotropic virus integration site 1 protein, EVI-1, MDS1 and EVI1 complex locus protein, Myelodysplasia syndrome 1 protein homolog, Mecom, Evi1, Mds1, Prdm3
Accession No.	P14404
Target Species	Mus musculus (Mouse)
Calculated MW	16.03 kDa
Target Length	Asp74-Ser193
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.