

Recombinant Mouse CASP8 Protein, N-His

Catalog No: #AP97914

Package Size: #AP97914-1 50ug #AP97914-2 100ug #AP97914-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Mouse CASP8 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Caspase-8, CASP-8, 3.4.22.61, Caspase-8 subunit p18, Caspase-8 subunit p10, Casp8
Accession No.	O89110
Target Species	Mus musculus (Mouse)
Calculated MW	30.04 kDa
Target Length	Gly236-Pro480
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.