

Recombinant Mouse GSTA3 Protein, N-His

Catalog No: #AP97915

Package Size: #AP97915-1 50ug #AP97915-2 100ug #AP97915-3 1mg

Description

Product Name	Recombinant Mouse GSTA3 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Glutathione S-transferase A3, 2.5.1.18, GST class-alpha member 3, Glutathione S-transferase Ya3, Glutathione S-transferase Yc, Gsta3, Gstyc
Accession No.	P30115
Target Species	Mus musculus (Mouse)
Calculated MW	27.52 kDa
Target Length	Met1-Ser221
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.