

Recombinant Mouse VTCN1/B7-H4 Protein, N-His

Catalog No: #AP97920

Package Size: #AP97920-1 50ug #AP97920-2 100ug #AP97920-3 1mg

Description

Product Name	Recombinant Mouse VTCN1/B7-H4 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	V-set domain containing T-cell activation inhibitor 1, B7 homolog 4, B7-H4, Immune costimulatory protein B7-H4, Protein B7S1, T cell costimulatory molecule B7x, Vtcn1, B7h4
Accession No.	Q7TSP5
Target Species	Mus musculus (Mouse)
Calculated MW	27.76 kDa
Target Length	Leu25-Asn255
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.