

Recombinant Mouse TICAM1 Protein, N-His

Catalog No: #AP97922

Package Size: #AP97922-1 50ug #AP97922-2 100ug #AP97922-3 1mg

Description

Product Name	Recombinant Mouse TICAM1 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	TIR domain-containing adapter molecule 1, TICAM-1, Toll-interleukin-1 receptor domain-containing adapter protein inducing interferon beta, TIR domain-containing adapter protein inducing IFN-beta, Ticam1, Trif
Accession No.	Q80UF7
Target Species	Mus musculus (Mouse)
Calculated MW	15.34 kDa
Target Length	Gln394-His510
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.