

Recombinant Mouse PHPT1 Protein, N-GST & C-His

Catalog No: #AP97929

Package Size: #AP97929-1 50ug #AP97929-2 100ug #AP97929-3 1mg

Description

Product Name	Recombinant Mouse PHPT1 Protein, N-GST & C-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	14 kDa phosphohistidine phosphatase, 3.9.1.3, Phosphohistidine phosphatase 1, PHPT1, Protein histidine phosphatase, PHP, Phpt1, Php14
Accession No.	Q9DAK9
Target Species	Mus musculus (Mouse)
Calculated MW	37.39 kDa
Target Length	Ala2-Asp84
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.