

Recombinant Mouse HDAC9 Protein, N-His

Catalog No: #AP97932

Package Size: #AP97932-1 50ug #AP97932-2 100ug #AP97932-3 1mg

Description

Product Name	Recombinant Mouse HDAC9 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Histone deacetylase 9, HD9, 3.5.1.98, Histone deacetylase 7B, HD7b, Histone deacetylase-related protein, MEF2-interacting transcription repressor MITR, Hdac9, Hdac7b, Hdrp, Mitr
Accession No.	Q99N13
Target Species	Mus musculus (Mouse)
Calculated MW	17.29 kDa
Target Length	Val37-Pro161
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.