

Recombinant Mouse ERRFI1 Protein, N-GST & C-His

Catalog No: #AP97935

Package Size: #AP97935-1 50ug #AP97935-2 100ug #AP97935-3 1mg

Description

Product Name	Recombinant Mouse ERRFI1 Protein, N-GST & C-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	ERBB receptor feedback inhibitor 1, Mitogen-inducible gene 6 protein homolog, MIG-6, Errf1, Mig6
Accession No.	Q99JZ7
Target Species	Mus musculus (Mouse)
Calculated MW	35.87 kDa
Target Length	Ser329-Arg399
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.