

Recombinant Mouse TAPBPL Protein, N-His

Catalog No: #AP97937

Package Size: #AP97937-1 50ug #AP97937-2 100ug #AP97937-3 1mg

Description

Product Name	Recombinant Mouse TAPBPL Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Tapasin-related protein, TAPASIN-R, TAP-binding protein-like, TAP-binding protein-related protein, TAPBP-R, Tapasin-like, Tapbpl
Accession No.	Q8VD31
Target Species	Mus musculus (Mouse)
Calculated MW	42.82 kDa
Target Length	Gln33-Pro406
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.