

Recombinant Mouse Ly6c2 Protein, N-GST & C-His

Catalog No: #AP97939

Package Size: #AP97939-1 50ug #AP97939-2 100ug #AP97939-3 1mg

Description

Product Name	Recombinant Mouse Ly6c2 Protein, N-GST & C-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Lymphocyte antigen 6C2, Ly-6C2, Ly6c2, Ly6c
Accession No.	P0CW03
Target Species	Mus musculus (Mouse)
Calculated MW	37.24 kDa
Target Length	Gly26-Ala108
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.