

Recombinant Rat CD13/ANPEP/APN Protein, N-His

Catalog No: #AP97941

Package Size: #AP97941-1 50ug #AP97941-2 100ug #AP97941-3 1mg

Description

Product Name	Recombinant Rat CD13/ANPEP/APN Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Aminopeptidase N, AP-N, pAPN, 3.4.11.2, Alanyl aminopeptidase, Aminopeptidase M, AP-M, Microsomal aminopeptidase, gp130, CD13, ANPEP
Accession No.	P15145
Target Species	Rattus norvegicus (Rat)
Calculated MW	25.37 kDa
Target Length	Thr75-Ser280
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.