

Recombinant ASFV p54/pE183L Protein, C-His

Catalog No: #AP97944

Package Size: #AP97944-1 50ug #AP97944-2 100ug #AP97944-3 1mg

Description

Product Name	Recombinant ASFV p54/pE183L Protein, C-His
Host Species	Mammalian Cells
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Inner membrane protein p54, pE183L, p54, E183L, E183L CDS, ASFV-Georgia_4-149
Accession No.	B2MUM6
Target Species	African swine fever virus (ASFV)
Calculated MW	16.87 kDa
Target Length	Ser53-Leu184
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.