

Recombinant ASFV P30/CP204L Protein, N-Fc

Catalog No: #AP97945

Package Size: #AP97945-1 50ug #AP97945-2 100ug #AP97945-3 1mg

Description

Product Name	Recombinant ASFV P30/CP204L Protein, N-Fc
Host Species	Mammalian Cells
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	p30, CP204L
Accession No.	AHN94305.1
Target Species	African swine fever virus (ASFV)
Calculated MW	52.61 kDa
Target Length	Met1-Ile201
Tag Info	Fc-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.