

Recombinant ASFV UBC/pI215L Protein, N-His

Catalog No: #AP97948

Package Size: #AP97948-1 50ug #AP97948-2 100ug #AP97948-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant ASFV UBC/pI215L Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Ubiquitin-conjugating enzyme E2, 2.3.2.23, E2 ubiquitin-conjugating enzyme, UBCv1, Ubiquitin carrier protein, Ubiquitin-protein ligase, pI215L, UBC, BA71V-144, I215L
Accession No.	P27949
Target Species	African swine fever virus (ASFV)
Calculated MW	25.01 kDa
Target Length	Met1-Thr200
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.